

ProBright Sonic Toothbrush - User Manual

Product Overview

The ProBright Sonic Toothbrush is a state-of-the-art dental hygiene tool designed to provide superior cleaning and comfort. Featuring advanced sonic technology, it offers 31,000 brush strokes per minute, ensuring efficient plaque removal and a gentle brushing experience. The ergonomic handle and intuitive control buttons make it both stylish and user-friendly.

Specifications

- **Model Number:** PBST1000
- **Dimensions:** 9.5 inches x 1.5 inches
- **Weight:** 0.5 lbs
- **Material:** BPA-free plastic
- **Battery:** Rechargeable Lithium-ion, providing up to 4 weeks of usage on a single charge
- **Brush Modes:** Clean, Sensitive, Whitening, Gum Care

Compatibility

The ProBright Sonic Toothbrush is compatible with the ProBright CS1000 and CS2000 charging stations. It also fits universal brush heads designed for sonic toothbrushes.

Usage Instructions

1. **Charging:** Place the toothbrush on the compatible charging station. Charge for at least 24 hours before first use.
2. **Attaching Brush Head:** Align the brush head with the handle and press down until it clicks into place.
3. **Turning On/Off:** Press the power button to turn on. Press again to switch between modes. Hold for 2 seconds to turn off.
4. **Brushing:** Apply a small amount of toothpaste, position the brush at a 45-degree angle, and gently glide across teeth.

5. **Timer:** Equipped with a 2-minute timer with 30-second interval reminders for even brushing.

Maintenance

- **Cleaning:** Rinse the brush head and handle under water. Avoid submersion.
- **Replacement:** Replace brush heads every 3 months or when bristles show wear.

Troubleshooting

- **Charging Issues:** Ensure the charging base is connected to power. Check contacts are clean.
- **Powering On/Off:** Avoid holding multiple buttons simultaneously; check if the battery needs recharging.
- **Mode Not Switching:** Reset the device by holding the power button for 10 seconds.

Customer Support

For assistance, contact ProBright Support: - **Phone:** 1-800-555-0199 - **Email:** support@probright.com - **Hours:** Monday-Friday, 9 AM - 5 PM EST - **Warranty:** 2-year limited warranty covering manufacturing defects.

FAQs

- **How often should I charge the toothbrush?** Typically, once every 3-4 weeks, depending on usage.
- **Is it safe to take the toothbrush on a plane?** Yes, it is airline-approved for both carry-on and checked luggage.
- **Can children use this toothbrush?** Suitable for children aged 7 and above under adult supervision.

By following this manual, you can ensure optimal performance and longevity of your ProBright Sonic Toothbrush.